


FNSGp(PGM) **1PM**; **2D**; **2c**; **3U**; **4G**; **5S**; **6Oct**; **7Rot**; **8Gp**; **9Ama**; **10Tm**; **11Azd**



O2k and DatLab file: P___ (A / B)		2017-		Operator:				
Sample type:		Cohort:		Sample code:				
Sample.Subsample number:		Unit:		Concentration:				
Medium: MiR06Cr								
Event	Mark name	State	Final conc. 2 ml O2k	Stock [mM]	Comment	Tit. [μl]	A	B
			~200 μM O2		MiR			
1ce	1ROX				Permeabilized cells			
1P			5 mM	2000		5		
1M	1PM	N(PM) _L	2 mM	400	Immediately after P	10		
2D	2D	N(PM) _P	1 / 2.5 mM	500		4 / 10		
2c	2c	N(PM) _{PC}	10 μM	4		5		
3U*	3U	N(PM) _E	Δ0.5 μM	1	CCCP	Δ1		
4G	4G	N(PGM) _E	10 mM	2000		10		
5S	5S	NS(PGM) _E	50 mM	1000		100		
6Oct	6Oct	FNS(PGM) _E	0.5 mM	100		10		
7Rot	7Rot	S _E	0.5 μM	1		1		
8Gp	8Gp	SGp _E	10 mM	1000		20		
9Ama	9Ama	ROX	2.5 μM	5		1		
10As			2 mM	800		5		
10Tm			0.5 mM	200	Immediately after As	5		
100					~20 min open			
10C	10Tm	Tm _E			Close and wait 5 min			
11Azd	11Azd	ROX	≥100 mM	4000	~5 min	100		